




1 aspirate supernatant, wash each well once with 200µl wash buffer




2 add 200 µl/well ELISA buffer, incubate for 1h




3 aspirate supernatant, wash each well once with 200µl wash buffer




4 transfer 100µl from each well of the assay plate to the ELISA plate, incubate for 2h




5 aspirate supernatant, wash each well once with 200µl wash buffer




6 add 100 µl/well detection anti-body, incubate for 1h




7 aspirate supernatant, wash each well once with 200µl wash buffer




8 add 100 µl/well Streptavidin-HRP, incubate for 30min



9 aspirate supernatant, wash each well thrice with 200µl wash



10 add 100 µl/well TBM solution, incubate for 15min in the dark



11 add 100 µl/well stop solution



12 measure OD at 450nm and 570nm for reference

day III: ELISA Assay


Perform each incubation step on a plate shaker.

- Prepare 150ml of wash buffer (0.05% Tween20 in PBS).
- Dilute ELISA Buffer 1:5 to 50ml in ddH₂O.
- Aspirate wells of the ELISA plate and wash once with 200 µl/well of wash buffer. 1
- Add 200µl ELISA buffer to each well of the ELISA plate. Incubate for 1h at room temperature. 2
- Aspirate wells and wash once with 200 µl/well of wash buffer. 3
- Transfer 100µl from the Assay plate to corresponding wells of the ELISA plate. Incubate for 2h at room temperature. 4
- Dilute the detection antibody 1:250 in ELISA buffer (42µl + 10.458ml).
- Aspirate the supernatant of the ELISA plate and wash once with 200µl washing buffer. Let the buffer soak in for 1min. 5
- Add 100µl detection antibody to each well and incubate for 1h at room temperature. 6
- Dilute Streptavidin-HRP 1:100 in ELISA buffer (110µl + 10.890ml).
- Aspirate the supernatant and wash once with 200µl washing buffer. Let the buffer soak in for 1min. 7
- Add 100µl Streptavidin-HRP to each well and incubate for 30min at room temperature. 8
- Aspirate the supernatant and wash thrice with 200µl washing buffer. Let the buffer soak in for 1min during every step. 9
- Add 100µl TMB solution to each well and incubate for 15min at room temperature while shaking in the dark. 10
- Add 100µl of stop solution (1M H₂SO₄) to each well. 11
- Measure OD at 450nm and 570nm for reference. 12


For analysis and acceptance criteria, please refer to the product information of the instaCELL pyrogen detection kit.

instaCELL® pyrogen detection kit (MAT) protocol







thaw cells for **3min** at **37°C**, dilute in **8ml** recovery buffer




dispense cells **100 µl/well**, incubate for **24h**




prepare coating solution



centrifuge for **3min** at **200xg**, resuspend in **10ml** assay buffer



add **100µl** of standard series and sample to each corresponding well, incubate for **18h**



add coating solution **100 µl/well**, incubate overnight at **4°C**

The instaCELL pyrogen detection kit can be conducted in three different versions:

- A. Quantitative Test
- B. Semi-Quantitative Test
- C. Reference Lot Comparison Test

This protocol refers to the Quantitative Test. For information on methods B, C and the product-specific validation to determine the maximum valid dilution (MVD), please refer to the product information on: www.accelerate.me/support/downloads


day I: preparation of cells

- Keep the cells on dry ice before thawing and process quickly.
- Add **60µl** of 50mM beta-mercaptoethanol to the assay buffer.
- Equilibrate the assay and recovery buffer to **37°C**.
- Thaw one vial of Assay Ready Cells in a water bath at **37°C** for **3min**. Prepare **8ml** of recovery buffer in a **15ml** centrifugation tube. Dispense the cells completely into the prepared tube. Rinse the cryovial once with **1ml** recovery buffer. ①
- Centrifuge for **3min** at **200xg** and carefully aspirate the supernatant. Resuspend the cell pellet in **10ml** of assay buffer. ②
- Dispense **100µl** of the cell suspension into each well of the provided assay plate. ③
- Incubate for **24h** in a humidified incubator at **37°C** and **5% CO₂**.

day II: preparation of test chemicals and ELISA plate

- Dilute the BRP stock solution **1:25** in assay buffer to obtain the BRP working solution. Dilute the working solution **1:50** to obtain **C1** of the BRP standard curve and perform 2-fold serial dilutions for a total of 7 dilutions (see plate layout).
- Each sample is tested in three concentrations. Use the concentration previously determined in the test for interfering factors as the highest sample concentration and prepare two-fold serial dilutions for concentrations 2 and 3, not exceeding the MVD.
- Prepare your sample in assay buffer at **2x** of the final concentration.
- Dilute the BRP working solution **1:4** to obtain the spiking solution. Add **10µl** of spiking solution to **990µl** of the test sample dilutions to spike your samples.
- Add **100µl** of BRP standard series, samples, spiked samples and blank in quadruplicates to the respective wells of the assay plate. Use Assay Buffer for the blank. Incubate for **18h** in a humidified incubator at **37°C** and **5% CO₂**. ④
- Dilute coating buffer **1:10** in cell culture water (**1.1ml + 9.9ml**). Dilute capture antibody **1:250** in coating buffer (**42µl + 10458µl**). ⑤
- Coat the provided 96-well ELISA plate by adding **100µl** of coating solution to each well. Incubate overnight at **4°C**. ⑥

kit content



not provided

DPBS
Tween 20
1M H₂SO₄
50mM β-mercaptoethanol

storage

Store Assay Ready Cells in liquid nitrogen (below -140°C)
Store all reagents and media at temperatures indicated on the label

limited product warranty

This warranty limits our liability to replace this product. acCELLerate does not provide any other warranties of any kind, expressed or implied. Without limitation, this includes implied warranties of merchantability or fitness for a particular purpose. acCELLerate shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product. The product is for research only and must not be used for diagnostic or therapeutic purposes.

limited use license

Assay Ready Cells are provided under a limited use license and shall only be used instantly for the particular purpose of the instaCELL bioassay kit. No further expansion, passaging or refreezing of the cells is permitted. When breaking the sealed bag of Assay Ready Cells, the user is explicitly accepting the terms of this limited use license. instaCELL® is a registered trademark of acCELLerate GmbH.

EUROPEAN OFFICE

+49 (160) 987 577 56

acCELLerate GmbH
Osterfeldstraße 12-14
22529 Hamburg
Germany

US OFFICE

+1 (732) 698 34 04

acCELLerate GmbH
400 Route 518
Skillman, NJ 08558
USA



please@accelerate.me
www.accelerate.me