


*ALWAYS READY,
ALWAYS **PERFORMING!***




RESEARCH

TRC – Transfection Ready Cells


- simple to use like a reagent
- transfect instantly after thawing
- work with all transfection protocols
- always available for immediate transfection
- reproducible vial after vial after vial




thaw



transfect



express



assay



HEALTHCARE



COSMETICS



FOOD



CHEMICALS



RESEARCH



SAFETY



QUALITY



DIAGNOSTICS

TRC – Transfection Ready Cells

mammalian cells ready for transfection

TRC - Transfection Ready Cells are cryopreserved cells which can be used for transfection **instantly after thawing**. No cultivation or passaging is required. TRCs work with basically all transfection protocols providing a **better reproducibility** due to a **more consistent** quality of the cells. Different transfection host cells are available which can be used for **cell based assays** or **transient protein production**.

TRC – Transfection Ready Cells

cat-no	cell line	unit
RE011-05	HEK293	5 vials á 1E7 cells
RE012-05	CHO-K1	5 vials á 1E7 cells
RE013-05	MDCK	5 vials á 1E7 cells
RE014-05	Vero	5 vials á 1E7 cells

If you are looking for any different cell types please call acCELLerate's customer service. We are constantly expanding of product line of Transfection Ready Cells.




Fig: TRC – Transfection Ready Cells were thawed at the day of transfection and washed once in buffer before they were instantly used for reverse transfection. Fresh cells from a continuously growing culture have been seeded the day before transfection. For transfection a luciferase expression plasmid was mixed with “Viromer Yellow” transfection reagent (Lipocalyx) and used according to the manufacturers advice. Luciferase expression was quantified 24 hours after transfection using Firefly Luc One-Step Glow Assay (Pierce).

EUROPEAN OFFICE

+49 (160) 98 75 77 56

acCELLerate GmbH
Osterfeldstraße 12-14
22529 Hamburg · Germany

please@accellerate.me
www.accellerate.me

US OFFICE

+1 (908) 262 25 92

acCELLerate GmbH
1 Jill Court, Building 16 Unit 10
Hillsborough, NJ 08844 · USA

