

Plate layout of the instaCELL® KeratinoSens® assay kit

	1	2	3	4	5	6	7	8	9	10	11	12
A	comp. 1 C12	comp. 1 C11	comp. 1 C10	comp. 1 C9	comp. 1 C8	comp. 1 C7	comp. 1 C6	comp. 1 C5	comp. 1 C4	comp. 1 C3	comp. 1 C2	comp. 1 C1
B	comp. 1 C12	comp. 1 C11	comp. 1 C10	comp. 1 C9	comp. 1 C8	comp. 1 C7	comp. 1 C6	comp. 1 C5	comp. 1 C4	comp. 1 C3	comp. 1 C2	comp. 1 C1
C	comp. 1 C12	comp. 1 C11	comp. 1 C10	comp. 1 C9	comp. 1 C8	comp. 1 C7	comp. 1 C6	comp. 1 C5	comp. 1 C4	comp. 1 C3	comp. 1 C2	comp. 1 C1
D	comp. 2 C12	comp. 2 C11	comp. 2 C10	comp. 2 C9	comp. 2 C8	comp. 2 C7	comp. 2 C6	comp. 2 C5	comp. 2 C4	comp. 2 C3	comp. 2 C2	comp. 2 C1
E	comp. 2 C12	comp. 2 C11	comp. 2 C10	comp. 2 C9	comp. 2 C8	comp. 2 C7	comp. 2 C6	comp. 2 C5	comp. 2 C4	comp. 2 C3	comp. 2 C2	comp. 2 C1
F	comp. 2 C12	comp. 2 C11	comp. 2 C10	comp. 2 C9	comp. 2 C8	comp. 2 C7	comp. 2 C6	comp. 2 C5	comp. 2 C4	comp. 2 C3	comp. 2 C2	comp. 2 C1
G	blank	PC C5	PC C4	PC C3	PC C2	PC C1	PC C5	PC C4	PC C3	PC C2	PC C1	blank
H	SC	SC	SC	SC	SC	SC	PC C5	PC C4	PC C3	PC C2	PC C1	blank

SC – solvent control

PC – positive Control

blank – no cells