


*READY
WHEN YOU ARE!*


frozen instant cells for bioassays

- + custom assay ready cells from sponsor cell lines
- + prepared for instant use, no cell expansion required
- + improved assay precision and reliability
- + lot size of up to 1,000 vials in custom aliquots
- + full documentation for GxP compliant assays
- + quick turnaround time



accuracy

the potency acquired with assay ready frozen instant cells very well correlates with results from continuously cultured cells. 45 different bioassays of various types have been compared.



precision

four independent batches of frozen instant cells from a recombinant target receptor cell line provide reproducible results in a GPCR signaling assay (IP-One htrf).



alignment of client's SOPs
with acCELLerate's process

optimization of cell expansion
and freezing protocol

pilot batch and definition
of product specifications

approval of production
process by client

manufacturing of
assay batch



large scale production



automated filling



controlled rate freezing



frozen instant cells

make your life easier


- + safe laboratory resources, no maintenance cell culture required
- + obtain reproducible results from reliably performing cells
- + flexibly schedule your assay, cells are instantly available at any time



thaw



dispense



add samples



incubate



measure

EUROPEAN OFFICE

+49 (160) 987 577 56

acCELLerate GmbH
Osterfeldstraße 12-14
22529 Hamburg · Germany

please@accellerate.me
www.accellerate.me

US OFFICE

+1 (732) 698 3404

acCELLerate GmbH
1 Jill Court, Bldg. 16 Unit 10
Hillsborough, NJ 08844 · USA

